

DSC-800

Differential scanning calorimeter



Differential scanning calorimetry (DSC), as a classical thermal analysis method to control the thermal effect at programmed temperature, has been widely used in research and development, process optimization, quality control and failure analysis in various fields of materials and chemistry. Using DSC method, we can study the phase transition of inorganic materials, the melting and crystallization process of polymer materials, the polymorphism of drugs, the solid/liquid ratio of foods such as oils and fats, etc.

main feature

- Simple operation, no testing experience, only a small amount of training.
- The software is suitable for computer screens with different resolutions.
- Double temperature probes ensure high accuracy and repeatability.
- During the experiment, there is no need for personnel supervision.
- Digital gas mass flowmeter automatically switches two gas flows.
- Software can be upgraded online for free.
- Seven-inch LCD display, the map and curve are clear at a glance.

technical parameter

| | | | |
|--------------|-------------|---------------------------|---|
| DSC range | 0~±200mW | data interface | USB standard interface and supporting software |
| heating rate | 0.1~50℃/min | temperature resolution | 0.01℃ |


| | | | |
|-------------------------------------|---|---------------------------|--|
| Temperature accuracy | ±0.1℃ | Temperature repeatability | ±0.1℃ |
| DSC accuracy | ±2% | DSC resolution | 0.001mW |
| DSC resolution | 0.001mW | Temperature control mode | Full automatic control |
| Curve scanning | Heating scanning | Atmosphere control | Gas mass flowmeter |
| display mode | 24bit color, 7-inch touch screen | working power supply | AC220V 50Hz/60Hz |
| Overall dimension (w*d*h) | 42*38*27cm | Net weight | 17KG |
| Packing size (w*d*h) | 58*45*40 | Gross weight | 20KG |
| temperature range | Room temperature ~ 800℃ | | |
| Computer configuration requirements | Minimum hardware configuration <ul style="list-style-type: none"> ● Celeron dual core ● 2GB system memory ● Resolution 1366*768 pixels or higher. ● Hard disk 500G | | Minimum software configuration <ul style="list-style-type: none"> ● Win10/11 ● Adobe PDF reader |

application example


Measure physical and chemical changes related to heat, such as glass transition temperature, melting point, melting temperature, crystallization and crystallization heat, phase transition reaction heat, thermal stability of products, curing/crosslinking, oxidation induction period, reaction kinetics, specific heat, etc. Note: the thermal stability test during oxidation induction period is applicable to G B/T 1 7 3 9 1-1 9 9 8.

Instrument principle


Differential scanning calorimetry (DSC) is a technique to measure the power difference between substance and reference under the control of programmed temperature. DSC series differential scanning calorimeter is mainly composed of heating furnace, host, microvolt amplifier, A/D converter, data acquisition system, gas flow control system, computer, printer and other components, supplemented by the switching of two channels of atmosphere, and the measurement results are processed by computer data processing system.




Experimental atlas




Interface of automatic calculation results of glass transition temperature



Enthalpy variable, epitaxial starting point, melting point calculation result interface



Enthalpy change ratio calculation result interface



Interface of calculation results of oxidation induction period

List of accessories for DSC-800

| serial number | Nominal name | quantity |
|---------------|---|----------|
| 1 | Heating furnace and main cabinet | one |
| 2 | U disk of differential scanning calorimeter working program | one |
| 3 | feed cable | one |
| 4 | signal cable | one |
| 5 | Glass tube fuse (3A) | four |
| 6 | tweezers | one |
| 7 | Sample spoon | one |



| | | |
|----|--|-----------|
| 8 | certificate | one |
| 9 | Pressure reducing valve connector | 2 |
| 10 | Ventilation joint | four |
| 11 | Aluminium crucible | 200 |
| 12 | Alumina crucible | 200 |
| 13 | Standard samples (In, Sn, Zn, Al,) | 0.5g each |
| 14 | Instructions for differential scanning calorimeter | one |
| 15 | Ventilation plastic pipe (blue) | 3 meters |
| 16 | Ventilation plastic pipe (orange) | 3 meters |
| 17 | Software encryption dog | one |