

YKGWF-D1 Particle Analyzer



Main Features:

- Meet the requirements of national standards of China Pharmacopeia、America Pharmacopeia、 Pharmaceutical package standard、 infusion devices GB8368.
- Full volume test, simulated clinical infusion process with objective and true data.
- Original packing and completely closed sampling to ensure no pollution in the process.
- It adopts new type of injection pump to ensure continuous sampling, high precision and fast speed.
- The liquid level automatic recognition function can distinguish the end point test without setting the sampling volume.
- The high performance laser light source and light energy compensation circuit, ensure the test precision of various samples of colorless, colored sample.
- The high performance processor, nearly 10000 counting channels, realize high precision data acquisition.
- The fully automatic calibration, built-in kinds of calibration curves, meet the requirements of different standards.
- The multi-level rights management for data security.
- Android operating system, built-in huge storage, data statistical analysis function.
- You can access the laboratory operation platform, network management.

Main Technical Index:

- Channel Set up:64 Channel, thousands of different particle sizes can be free set up, accuracy 0.1 μ m
- Test Range:1~500 μ m



0086 16601757347
inquiry@yukelab.com
www.yukelab.com
0086 021 59570209

- Count Range:0~9999999 particle
- Sample Volume:1ml~full volume, accuracy 0.1ml
- Accuracy of Sample Volume: $< \pm 0.5\%$
- Sample Speed: 5~80ml/min, arbitrary set up
- Counting accuracy: $< \text{required value} \pm 5\%$
- Channel Resolution: $>95\%$
- Relative Standard Deviation: $\text{RSD} < 1.5\%$ (standard particle ≥ 1000 particle/ml)
- The best detection concentration: 0~18000 particle/ml
- Mixing Speed: 0~1000r/min, impeller mixer of adjustable
- Work Temperature:10~40°C
- Power: AC220V $\pm 10\%$; 50Hz; $<80\text{W}$
- Data output: Built-in Pin mini-printer; RS232 interface can connect computer; USB saved interface.